

Product Datasheet

Human EDIL3 protein (orb391640)

Description

Recombinant human EDIL3 protein is

system and the target gene e...

Reactivity

Human

Endotoxins

Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test.

Tested

ELISA, MS, SDS-PAGE, WB

Applications
Target

EDIL3

Preservatives

Liquid in a 0.2 μM filtered solution of 20mM Tris-HCl, 300mM NaCl, pH 8.0.

Form/Appearance

Liquid in a 0.2 μM filtered solution of 20mM Tris-HCl, 300mM NaCl, pH 8.0.

Storage

Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 μg/ml. Dissolve the lyophilized protein in 1X PBS. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

Protein Sequence

DIC DPNPCENGGI CLPGLADGSF SCECPDGFTD
PNCSSVVEVA SDEEPTSAG PCTPNPCHNG
GTCEISEAYR GDTFIGYVCK CPRGFNGIHC
QHNINECEVE PCKNGGICTD LVANYSCECP
GEFMGRNCQY KCSGPLGIEG GIISNQQITA
SSTHRALFGL QKWYPYARL NKKGLINAWT
AAENDRWPWI QINLQRKMRV TGVITQGAKR
IGSPEYIKSY KIAYSNDGKT WAMYKVKGTN
EDMVFRGNID NNTPYANSFT PPIKAQYVRL
YPQVCRRHCT LRMELLGCEL SGCSEPLGMK
SGHIQDYQIT ASSIFRTLNM DMFTWEPRKA
RLDKQGKVNA WTSGHNDQSQ WLQVDLLVPT
KVTGIITQGA KDFGHVQFVG SYKLAYSNDG
EHWTVYQDEK QRKDKVFQGN FDNDTHRKNV
IDPPIYARHI RILPWSWYGR ITLRSELLGC
TEEEVDHHHH HH

Purity

Greater than 95% as determined by SEC-HPLC and reducing SDS-PAGE.