

Product Datasheet

Human EFNB1 protein (orb391649)

Description

Recombinant human EFNB1 protein is

system and the target gene e...

Reactivity

Human

Endotoxins

Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test.

Tested

ELISA, MS, SDS-PAGE, WB

Applications

Target

EFNB1

Preservatives

Lyophilized from a 0.2 μM filtered solution of 20mM PB, 150mM NaCl, pH 7.4.

Form/Appearance

Lyophilized from a 0.2 μM filtered solution of 20mM PB, 150mM NaCl, pH 7.4.

Storage

Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 2-8°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.

Note

For research use only

Application notes

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 μg/ml. Dissolve the lyophilized protein in 1X PBS. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

Protein Sequence

LAKNLEPVSW SSLNPKFLSG KGLVIYPKIG
DKLDIICPRA EAGRPYEYYK LYLVRPEQAA
ACSTVLDPNV LVTCLRPEQE IRFTIKFQEF
SPNYMGLEFK KHHYYITST SNGSLEGLN
REGGVCRTTR MKIIMKVGQD PNAVTPEQLT
TSRPSKEADN TVKMATQAPG SRGSLGSDSG
KHETVNQEEK SGPGASGGSS GDPDGVDHHH
HHH

Purity

Greater than 95% as determined by reducing SDS-PAGE.

Source

Human cells

Uniprot ID

[P98172](#)

Entrez

[1947](#)

Expiration Date

6 months from date of receipt.