

Product Datasheet

Human F11R protein (orb391692)

Description

Recombinant human F11R protein is

system and the target gene en...

Reactivity

Human

Endotoxins

Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test.

Tested

ELISA, MS, SDS-PAGE, WB

Applications
Target

F11R

Preservatives

Lyophilized from a 0.2 μM filtered solution of 20mM Tris-HCl, 50mM NaCl, 100mM Glycine, pH 7.5.

Form/Appearance

Lyophilized from a 0.2 μM filtered solution of 20mM Tris-HCl, 50mM NaCl, 100mM Glycine, pH 7.5.

Storage

Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 2-8°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.

Note

For research use only

Application notes

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 μg/ml. Dissolve the lyophilized protein in 1X PBS. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

Protein Sequence

SVTVHSSEPE VRIPENNPVK LSCAYSGFSS
PRVEWKFDQG DTTRLVCYNN KITASYEDRV
TFLPTGITFK SVTREDTGTY TCMVSEEGN
SYGEVKVKLI VLVPPSKPTV NIPSSATIGN
RAVLTCEQD GSPPSEYTWf KDGIVMPTNP
KSTRAFSNSS YVLNPTTGEL VFDPLSASDT
GEYSCEARNg YGTPMTSNAV RMEAVERNVG
VDHHHHHH

Purity

Greater than 95% as determined by SEC-HPLC and reducing SDS-PAGE.

Source

Human cells

Uniprot ID
[Q9Y624](#)
Entrez
[50848](#)