

Product Datasheet

Human FETUB protein (orb391719)

Description

Recombinant human FETUB protein is

system and the target gene e...

Reactivity

Human

Endotoxins

Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test.

Tested

ELISA, MS, SDS-PAGE, WB

Applications

Target

FETUB

Preservatives

Lyophilized from a 0.2 μM filtered solution of 20mM PB, 150mM NaCl, pH 7.2.

Form/Appearance

Lyophilized from a 0.2 μM filtered solution of 20mM PB, 150mM NaCl, pH 7.2.

Storage

Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 2-8°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.

Note

For research use only

Application notes

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 μg/ml. Dissolve the lyophilized protein in 1X PBS. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

Protein Sequence

CGAMSPPQLA LNPSALLSRG CNDSDVLAVA
GFALRDINKD RKDGYVLRN RVNDAQEYRR
GGLGSLFYLT LDVLETDCV LRKKAQDCG
MRIFESVYG QCKAIFYMNN PSRVLYLAAY
NCTLRPVSKK KIYMTCPDCP SSIPTDSSNH
QVLEAATESL AKYNNENTSK QYSLFKVTRA
SSQWVVGPSY FVEYLIKESV CTKSQASSCS
LQSSDSVPVG LCKGSLTRTH WEKFSVSTCD
FFESQAPATG SENSANQKP TNLPKVEESQ
QKNTPTDPS SKAGPRGSVQ YLPDLDDKNS
QKGPQEAFF VHLDLTTNPQ GETLDISFLF
LEPMEEKLVV LPFPKEKART AECGPAQNA
SPLVLPPVDH HHHHH

Purity

Greater than 95% as determined by reducing SDS-PAGE.

Source

Human cells