

Product Datasheet

Human FKBP7 protein (orb391749)

Description

Recombinant human FKBP7 protein is

system and the target gene e...

Reactivity

Human

Endotoxins

Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test.

Tested Applications

ELISA, MS, SDS-PAGE, WB

Target

FKBP7

Preservatives

Liquid in a 0.2 μM filtered solution of 20mM Tris-HCl, 150mM NaCl, 1mM CaCl₂, 10% Glycerol, pH 7.5.

Form/Appearance

Liquid in a 0.2 μM filtered solution of 20mM Tris-HCl, 150mM NaCl, 1mM CaCl₂, 10% Glycerol, pH 7.5.

Storage

Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 μg/ml. Dissolve the lyophilized protein in 1X PBS. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

Protein Sequence

QRQKKEESTE EVKIEVLHRP ENCSKTSKKG
DLLNAHYDGY LAKDGSKFYC SRTQNEGHPK
WFLVGVGQVI KGLDIAMTDM CPGEKRKQVI
PPSFAYGKEG YAEGKIPDA TLIFEIELYA
VTGPRSIET FKQIDMDNDR QLSKAEINLY
LQREFEKDEK PRDKSYQDAV LEDIFKKNDH
DGDGFISPKE YNVYQHDELV DHHHHHH

Purity

Greater than 95% as determined by reducing SDS-PAGE.

Source

Human cells

Uniprot ID

[Q9Y680](#)

Entrez

[51661](#)

Expiration Date

6 months from date of receipt.