

Product Datasheet

Human FLRT3 protein (orb391752)

Description

Recombinant human FLRT3 protein is

system and the target gene e...

Reactivity

Human

Endotoxins

Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test.

Tested

ELISA, MS, SDS-PAGE, WB

Applications

Target

FLRT3

Preservatives

Lyophilized from a 0.2 μM filtered solution of PBS, pH 7.4.

Form/Appearance

Lyophilized from a 0.2 μM filtered solution of PBS, pH 7.4.

Storage

Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 2-8°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.

Note

For research use only

Application notes

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 μg/ml. Dissolve the lyophilized protein in 1X PBS. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

Protein Sequence

KSCPSVCRCD AGFIYCND RF LTSIPTGIPE
DATTLYLQNN QINNAGIPSD LKNLLKVERI
YLYHNSLDEF PTNLPKYVKE LHLQENNIRT
ITYDSLSKIP YLEELHDDN SVSAVSIEEG
AFRDSNYLRL LFLSRNHLST IPWGLPRTIE
ELRLDDNRIS TISSPSLQGL TSLKRLVLDG
NLLNNHGLGD KVFFNLVNL ELSLVRNSLT
AAPVNLPGTN LRKLYLQDNH INRVPPNAFS
YLRQLYRLDM SNNNLSNLPQ GIFDDLNDIT
QLILRNNPWY CGCKMKWVRD WLQSLPVKVN
VRGLMCQAPE KVRGMAIKDL NAELFDCKDS
GIVSTIQITT AIPNTVYPAQ QWPAPVTKQ
PDIKNPKLTK DHQTTGSPSR KTITITVKS
TSDTIHISWK LALPMTALRL SWLKLGHSPA
FGSITETIVT GERSEYLTA LEPDSPYKVC
MVPMETSNLY LFDETPVCIE TETAPLRMYN
PTTTLNREQE KEPYKNPNLP VDHSHHHH

Biorbyt Ltd.

7 Signet Court, Swann's Road, Cambridge, CB5 8LA, United Kingdom

Email: info@biorbyt.com | Phone: +44 (0) 1223 859-353 | Fax: +44 (0)1223 280 240

Biorbyt LLC.

68 TW Alexander Drive
Research Triangle Park
Durham, North Carolina
27709, United States

Email: info@biorbyt.com | Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558