

Product Datasheet

Human FOLR1 protein (orb391757)

Description

Recombinant human FOLR1 protein is

system and the target gene e...

Reactivity

Human

Endotoxins

Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test.

Tested

ELISA, MS, SDS-PAGE, WB

Applications

Target

FOLR1

Preservatives

Lyophilized from a 0.2 μM filtered solution of PBS, pH 7.4.

Form/Appearance

Lyophilized from a 0.2 μM filtered solution of PBS, pH 7.4.

Storage

Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 2-8°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.

Note

For research use only

Application notes

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 μg/ml. Dissolve the lyophilized protein in 1X PBS. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

Protein Sequence

RIAWARTELL NVCMNAKHHK EKPGPEDKLH
EQCRPWKRKNA CCSTNTSQEA HKDVSYLRYF
NWNHCGEMAP ACKRHFIQDT CLYECSPNLG
PWIQQVDQSW RKERVNLNPL CKEDCEQWWE
DCRTSYTCKS NWHKGWNWTS GFNKA VGAA
CQPFHFYFPT PTVLCNEIWT HSYKVSNSR
GSGRCIQMWF DPAQGNPNEE VARFYAAAMS
VDHHHHHH

Purity

Greater than 95% as determined by reducing SDS-PAGE.

Source

Human cells

Uniprot ID

[P15328](#)

Entrez

[2348](#)

Expiration Date

6 months from date of receipt.