

Product Datasheet

Human GALM protein (orb391776)

Description

Recombinant human GALM protein is

and the target gene enco...

Reactivity

Human

Endotoxins

Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test.

Tested

ELISA, MS, SDS-PAGE, WB

Applications

Target

GALM

Preservatives

Liquid in a 0.2 μM filtered solution of 20mM PB, 150mM NaCl, pH 7.2.

Form/Appearance

Liquid in a 0.2 μM filtered solution of 20mM PB, 150mM NaCl, pH 7.2.

Storage

Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 μg/ml. Dissolve the lyophilized protein in 1X PBS. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

Protein Sequence

ASVTRAVFGE LPSGGGTVEK FQLQSDLLRV
DIISWGCTIT ALEVKDRQGR ASDVVLGFAE
LEGYLQKQPY FGAVIGRVAN RIAKGTFKVD
GKEYHLAINK EPNSLHGGVR GFDKVLWTPR
VLSNGVQFSR ISPDGEEGYP GELKVWVTYT
LDGGELIVNY RAQASQATPV NLTNHSYFNL
AGQASPNIND HEVTIEADTY LPVDETLIPT
GEVAPVQGT AFDLRKPVELG KHLQDFHLNG
FDHNFCLKGS KEKHFCARVH HAASGRVLEV
YTTQPGVQFY TGNFLDGTGK GKNQAVYPKH
SGFCLETQNW PDAVNQPRFP PVLLRPGEY
DHTTWFKFSV ALEHHHHHHH

Purity

Greater than 95% as determined by reducing SDS-PAGE.

Source

E. coli

Uniprot ID

[Q96C23](#)

Entrez

[130589](#)

Biorbyt Ltd.

7 Signet Court, Swann's Road, Cambridge, CB5 8LA, United Kingdom

Email: info@biorbyt.com | Phone: +44 (0) 1223 859-353 | Fax: +44 (0)1223 280 240

Biorbyt LLC.

68 TW Alexander Drive
Research Triangle Park
Durham, North Carolina
27709. United States

Email: info@biorbyt.com | Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558