

Product Datasheet

Human GTF2B protein (orb391827)

Description

Recombinant human GTF2B protein is

and the target gene enc...

Reactivity

Human

Endotoxins

Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test.

Tested

Applications

ELISA, MS, SDS-PAGE, WB

Target

GTF2B

Preservatives

Lyophilized from a 0.2 μM filtered solution of 20mM Tris-HCl, 100mM NaCl, pH 8.0.

Form/Appearance

Lyophilized from a 0.2 μM filtered solution of 20mM Tris-HCl, 100mM NaCl, pH 8.0.

Storage

Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 2-8°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.

Note

For research use only

Application notes

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 μg/ml. Dissolve the lyophilized protein in 1X PBS. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

Protein Sequence

MSPILGYWKI KGLVQPTRL LLEYLEEKYEE
HLYERDEGDK WRNKKFELGL EFPNLPYYID
GDVKLTQSMA IIRYADKHN MLGGCPKERA
EISMLEGAVL DIRYGVSRIA YSKDFETLKV
DFLSKLPEML KMFEDRLCHK TYLNGDHVTH
PDFMLYDALD VVLYMDPMCL DAFPKLVCFK
KRIEAIQID KYLKSSKYIA WPLQGWQATF
GGGDHPPKSD LVPRGSPEFH MASTSRLDAL
PRVTCPNHPD AILVEDYRAG DMICPECGLV
VGDRVIDVGS EWRTFSNDKA TKDPSRVGDS
QNPLSDGDL STMIGKGTGA ASFDEFGNSK
YQNRRTMSSS DRAMNFAKE ITTMADRINL
PRNIVDRNTN LFKQVYEKLS LKGRANDAIA
SACLYIACRQ EGVPRTFKEI CAVSRISKKE
IGRCFKLILK ALETSDVLIT TGD FMSR FCS
NLCLPKQVQM AATHIARKAV ELDLVPGRSP
ISVAAAAIYM ASQASAEKRT QKEIGDIAGV
ADVTIROS YR I IYPRAPDLE PTDFKEDTPV

Biorbyt Ltd.

7 Signet Court, Swann's Road, Cambridge, CB5 8LA, United Kingdom

Email: info@biorbyt.com | Phone: +44 (0) 1223 859-353 | Fax: +44 (0)1223 280 240

Biorbyt LLC.

68 TW Alexander Drive
Research Triangle Park
Durham, North Carolina
27709, United States

Email: info@biorbyt.com | Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558