

Product Datasheet

Human HEPACAM protein (orb391844)

Description

Recombinant human HEPACAM protein is

system and the target gene...

Reactivity

Human

Endotoxins

Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test.

Tested

ELISA, MS, SDS-PAGE, WB

Applications

Target

HEPACAM

Preservatives

Lyophilized from a 0.2 μM filtered solution of 20mM PB, 150mM NaCl, pH 7.2.

Form/Appearance

Lyophilized from a 0.2 μM filtered solution of 20mM PB, 150mM NaCl, pH 7.2.

Storage

Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 2-8°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.

Note

For research use only

Application notes

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 μg/ml. Dissolve the lyophilized protein in 1X PBS. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

Protein Sequence

VNITSPVRLI HGTVGKSALL SVQYSSTSSD
RPVVKWQLKR DKPVTVVQSI GTEVIGTLRP
DYRDRIRLFE NGSLLLSDLQ LADEGTYEVE
ISITDDTFTG EKTINLTVDV PISRPQVLVA
STTVLELSEA FTLNCSHENG TKPSYTWLKD
GKPLLNSRM LLSPDQKVLV ITRVLMEDDD
LYSCMVENPI SQGRSLPVKI TVYRRSSVDH
HHHHH

Purity

Greater than 95% as determined by SEC-HPLC and reducing SDS-PAGE.

Source

Human cells

Uniprot ID

[Q14CZ8](#)

Entrez

[220296](#)

Expiration Date

6 months from date of receipt.