

Product Datasheet

Human HTN3 protein (orb391869)

Description

Recombinant human HTN3 protein is

and the target gene enco...

Reactivity

Human

Endotoxins

Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test.

Tested

ELISA, MS, SDS-PAGE, WB

Applications
Target

HTN3

Preservatives

Liquid in a 0.2 μM filtered solution of PBS, pH 7.4.

Form/Appearance

Liquid in a 0.2 μM filtered solution of PBS, pH 7.4.

Storage

Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 μg/ml. Dissolve the lyophilized protein in 1X PBS. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

Protein Sequence

MSPILGYWKI KGLVQPTRL L LEYLEEKYEE
HLYERDEGDK WRNKKFELGL EFPNLPYYID
GDVKLTQSMA IIRYIADKHN MLGGCPKERA
EISMLEGAVL DIRYGVSRIA YSKDFETLKV
DFLSKLPEML KMFEDRLCHK TYLNGDHVTH
PDFMLYDALD VVLYMDPMCL DAFPKLVCFK
KRIEAIQID KYLKSSKYIA WPLQGWQATF
GGGDHPPKSD LVPRGSENL FQGHMDSHAK
RHHGYKRKFH EKHSHRGYR SNLYDN

Purity

Greater than 95% as determined by reducing SDS-PAGE.

Source

E. coli

Uniprot ID

P15516

Entrez

3347

Expiration Date

6 months from date of receipt.