

Product Datasheet

Mouse IFNGR1 protein (orb391894)

Description

Recombinant mouse IFNGR1 protein is

system and the target gene ...

Reactivity

Mouse

Endotoxins

Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test.

Tested Applications

ELISA, MS, SDS-PAGE, WB

Target

IFNGR1

Preservatives

Lyophilized from a 0.2 μM filtered solution of 20mM PB, 150mM NaCl, pH 7.4.

Form/Appearance

Lyophilized from a 0.2 μM filtered solution of 20mM PB, 150mM NaCl, pH 7.4.

Storage

Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 2-8°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.

Note

For research use only

Application notes

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 μg/ml. Dissolve the lyophilized protein in 1X PBS. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

Protein Sequence

ALTSTEDPEP PSVPVPTNVL IKSYNLNPVV
CWEYQNMSQT PFTVQVKVY SGSWTDSCN
ISDHCCNIYE QIMYPDVSAW ARVKAKVGQK
ESDYARSKEF LMCLKGKVG PGLIIRRKKE
EQLSVLVFHP EVVNGESQG TMFGDGSTCY
TFDYTVYVEH NRSGEILHTK HTVEKEECNE
TLCELNISVS TLDSRYCISV DGISSFQVQR
TEKSKDVCIP PFHDDRKDVD HHHHHH

Purity

Greater than 95% as determined by reducing SDS-PAGE.

Source

Human cells

Uniprot ID
P15261
Entrez
15979
Expiration Date

6 months from date of receipt.