

Product Datasheet

Human IL18RAP protein (orb391942)

Description

Recombinant human IL18RAP protein is

system and the target gene...

Reactivity

Human

Endotoxins

Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test.

Tested

ELISA, MS, SDS-PAGE, WB

Applications

Target

IL18RAP

Preservatives

Lyophilized from a 0.2 μM filtered solution of 20mM PB, 150mM NaCl, pH 7.2.

Form/Appearance

Lyophilized from a 0.2 μM filtered solution of 20mM PB, 150mM NaCl, pH 7.2.

Storage

Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 2-8°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.

Note

For research use only

Application notes

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 μg/ml. Dissolve the lyophilized protein in 1X PBS. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

Protein Sequence

FNISGCSTKK LLWTYSTRSE EEVLFCDLP
EPQKSHFCHR NRLSPKQVPE HLPFMGSNDL
SDVQWYQQPS NGDPLEDIRK SYPHIIQDKC
TLHFLTPGVN NSGSYICRPK MIKSPYDVAC
CVKMILEVKP QTNASCEYSA SHKQDLLLS
TGSISCPSLS CQSDAQSPAV TWYKNGKLLS
VERSNRIVVD EYDYHQGT YCDYTQSDTV
SSWTVRAVVQ VRTIVGDTKL KPDILDPVED
TLEVELGKPL TISCKARFGF ERVFNPIKW
YIKDSLEWE VSVPEAKSIK STLKDEIER
NIILEKVTQR DLRRKFVCFV QNSIGNTTQS
VOLKEKRVDD IEGRMDEPKS CDKTHTCPPC
PAPELLGGPS VFLFPPKPKD TLMISRTPEV
TCVVVDVSHE DPEVKFNWYV DGVEVHNAKT
KPREEQYNST YRVVSVLTVL HQDWLNGKEY
KCKVSNKALP APIEKTISKA KGQPREPQVY
TLPPSREEMT KNQVSLTCLV KGFYPSDIAV
FWESNGOPEN NYKTTTPVLD SDGSFEIYSK

Biorbyt Ltd.

7 Signet Court, Swann's Road, Cambridge, CB5 8LA, United Kingdom

Email: info@biorbyt.com | Phone: +44 (0) 1223 859-353 | Fax: +44 (0)1223 280 240

Biorbyt LLC.

68 TW Alexander Drive
Research Triangle Park
Durham, North Carolina
27709. United States

Email: info@biorbyt.com | Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558

Source

Human cells