

Product Datasheet

Human IL6 protein (orb391992)

Description

Recombinant human IL6 protein is produced

gene encodi...

Reactivity

Human

Endotoxins

Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test.

Tested

ELISA, MS, SDS-PAGE, WB

Applications

Target

IL6

Preservatives

Lyophilized from a 0.2 μM filtered solution of 20mM HAc-NaAc, 150mM NaCl, pH 5.0

Form/Appearance

Lyophilized from a 0.2 μM filtered solution of 20mM HAc-NaAc, 150mM NaCl, pH 5.0

Storage

Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 2-8°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.

Note

For research use only

Application notes

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 μg/ml. Dissolve the lyophilized protein in 1X PBS. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

Protein Sequence

PVPPGEDSKD VAAPHRQPLT SSERIDKQIR
YILDGISALR KETCNKSNMC ESSKEALAEN
NLNLPKMAEK DGCFQSGFNE ETCLVKIITG
LLEFEVYLEY LQNRFESSEE QARAVQMSTK
VLIQFLQKKA KNLDAITTPD PTTNASLLTK
LQAQNQWLQD MTTHLILRSF KEFLQSSLRA
LRQM

Purity

Greater than 98% as determined by RP-HPLC and reducing SDS-PAGE.

Source

E. coli

Uniprot ID

[P05231](#)

Entrez

[3569](#)

Expiration Date

6 months from date of receipt.