

Product Datasheet

AKAP3 antibody (orb39201)

Description	Peptide affinity purified rabbit polyclonal antibody.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	AKAP3
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	94751
Uniprot ID	O75969
NCBI	NP_001265238.1 , NP_006413.3
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

