

Product Datasheet

Human LRG1 protein (orb392090)

Description

Recombinant human LRG1 protein is

system and the target gene en...

Reactivity

Human

Endotoxins

Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test.

Tested

ELISA, MS, SDS-PAGE, WB

Applications

Target

LRG1

Preservatives

Liquid in a 0.2 μM filtered solution of 20mM Tris, 200mM NaCl, 10% Glycerol, pH 8.0.

Form/Appearance

Liquid in a 0.2 μM filtered solution of 20mM Tris, 200mM NaCl, 10% Glycerol, pH 8.0.

Storage

Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 μg/ml. Dissolve the lyophilized protein in 1X PBS. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

Protein Sequence

VTLSPKDCQV FRSDHGSSIS CQPPAEIPGY
LPADTVHLAV EFFNLTHLPA NLLQGASKLQ
ELHLSSNGLE SLSPEFLRPV PQLRVLDLTR
NALTGLPPGL FQASATLDTL VLKENQLEVL
EVSWLHGLKA LGHLDLSGSR LRLPPGLLA
NFTLLRTLDL GENQLETLPD DLLRGPLQLE
RLHLEGKQLQ VLGKDLLLPQ PDLRYLFLNG
NKLARVAAGA FQGLRQLDML DLSNNSLASV
PEGLWASLGQ PNWDMRDGFD ISGNPWICDQ
NLSDLRWLQ AQKDKMFSQN DTRCAGPEAV
KGQTLAVAK SQVDDIEGRM DEPKSCDKTH
TCPPCPAPEL LGGPSVFLFP PKPKDTLMIS
RTPEVTCVVV DVSHEDPEVK FNWYVDGVEV
HNAKTKPREE QYNSTYRVVS VLTVLHQDWL
NGKEYKCKVS NKALPAPIEK TISKAKGQPR
EPQVYTLPPS REEMTKNQVS LTCLVKGFYP
SDIAVEWESN GPENNYKTT PPVLDSDGSF
FLYSKLTVDK SRWQQGNVFS CSVMHEALHN
HYTQKSLSLS PGKHHHHHHH

Biorbyt Ltd.

7 Signet Court, Swann's Road, Cambridge, CB5 8LA, United Kingdom

Email: info@biorbyt.com | Phone: +44 (0) 1223 859-353 | Fax: +44 (0)1223 280 240

Biorbyt LLC.

68 TW Alexander Drive
Research Triangle Park
Durham, North Carolina
27709, United States

Email: info@biorbyt.com | Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558