

Product Datasheet

Human LYG2 protein (orb392108)

Description

Recombinant human LYG2 protein is

system and the target gene en...

Reactivity

Human

Endotoxins

Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test.

Tested

ELISA, MS, SDS-PAGE, WB

Applications

Target

LYG2

Preservatives

Liquid in a 0.2 μM filtered solution of 20mM PB, 150mM NaCl, 10% Glycerol, pH 7.2.

Form/Appearance

Liquid in a 0.2 μM filtered solution of 20mM PB, 150mM NaCl, 10% Glycerol, pH 7.2.

Storage

Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 μg/ml. Dissolve the lyophilized protein in 1X PBS. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

Protein Sequence

SYPFSSHMKP HLHPRLYHGC YGDIMTMKTS
GATCDANSVM NCGIRGSEMF AEMDLRAIKP
YQTLIKEVGQ RHCVDPAVIA AIISRESHGG
SVLQDGWDHR GLKFGLMQLD KQTYHPVGAW
DSKEHLSQAT GILTERIKAI QKKFPTWSVA
QHLKGGLSAF KSGIEAIATP SDIDNDFVND
IIARAKFYKR QSFVDHHHHH H

Purity

Greater than 95% as determined by reducing SDS-PAGE.

Source

Human cells

Uniprot ID

[Q86SG7](#)

Entrez

[254773](#)

Expiration Date

6 months from date of receipt.