

## Product Datasheet

### Human NOL3 protein (orb392196)

**Description**

Recombinant human NOL3 protein is

and the target gene enco...

**Reactivity**

Human

**Endotoxins**

Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test.

**Tested**

ELISA, MS, SDS-PAGE, WB

**Applications**
**Target**

NOL3

**Preservatives**

Liquid in a 0.2 μM filtered solution of 25mM Tris-HCl, 1mM DTT, 1mM EDTA, 2mM β-ME, 20% Glycerol, pH 7.5.

**Form/Appearance**

Liquid in a 0.2 μM filtered solution of 25mM Tris-HCl, 1mM DTT, 1mM EDTA, 2mM β-ME, 20% Glycerol, pH 7.5.

**Storage**

Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

**Note**

For research use only

**Application notes**

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 μg/ml. Dissolve the lyophilized protein in 1X PBS. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

**Protein Sequence**

MGNAQERPSE TIDRERKRLV ETLQADSGLL  
LDALLARGVL TGPEYEALDA LPDAERRVRR  
LLLLVQKGGE AACQELLRCA QRTAGAPDPA  
WDWQHVGPGY RDRSYDPPCP GHWTPEAPGS  
GTTCPLPRA SDPDEAGGPE GSEAVQSGTP  
EEPEPELEAE ASKEAPEPE PEPELEPEAE  
AEPEPELEPE PDPEPEPDFE ERDESEDS

**Purity**

Greater than 95% as determined by SEC-HPLC and reducing SDS-PAGE.

**Source**

E. coli

**Uniprot ID**

**060936**

**Entrez**

**8996**

**Expiration Date**

6 months from date of receipt.