

Product Datasheet

Human NOL3 protein (orb392197)

Description

Recombinant human NOL3 protein is

the target gene encod...

Reactivity

Human

Endotoxins

Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test.

Tested

ELISA, MS, SDS-PAGE, WB

Applications
Target

NOL3

Preservatives

Lyophilized from a 0.2 μM filtered solution of 20mM Tris-HCl, 1mM DTT, 150mM NaCl, pH 8.0.

Form/Appearance

Lyophilized from a 0.2 μM filtered solution of 20mM Tris-HCl, 1mM DTT, 150mM NaCl, pH 8.0.

Storage

Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 2-8°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.

Note

For research use only

Application notes

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 μg/ml. Dissolve the lyophilized protein in 1X PBS. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

Protein Sequence

MSPILGYWKI KGLVQPTRLL LEYLEEKYEE
HLYERDEGDK WRNKKFELGL EFPNLPYYID
GDVKLTQSMA IIRYADKHN MLGGCPKERA
EISMLEGAVL DIRYGVSRIA YSKDFETLKV
DFLSKLPEML KMFEDRLCHK TYLNGDHVTH
PDFMLYDALD VVLYMDPMCL DAFPKLVCFK
KRIEAIQID KYLKSSKYIA WPLQGWQATF
GGGDHPPKSD LVPRGSMGNA QERPSETIDR
ERKRLVETLQ ADSGLLLDAL LARGVLTGPE
YEALDALPDA ERRVRRLLLL VQGKGAAACQ
ELLCAQRTA GAPDPAWDWQ HVGPGYRDRS
YDPPCPGHWT PEAPGSGTTC PGLPRASDPD
EAGGPEGSEA VQSGTPEEPE PELEAEASKE
AEPEPEPEPE LEPEAEAEPE PELEPEPDPE
PEPDFEERDE SEDS

Biorbyt Ltd.

7 Signet Court, Swann's Road, Cambridge, CB5 8LA, United Kingdom

Email: info@biorbyt.com | Phone: +44 (0) 1223 859-353 | Fax: +44 (0)1223 280 240

Biorbyt LLC.

68 TW Alexander Drive
Research Triangle Park
Durham, North Carolina
27709, United States

Email: info@biorbyt.com | Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558