

Product Datasheet

Human NTNG1 protein (orb392214)

Description

Recombinant human NTNG1 protein is

system and the target gene e...

Reactivity

Human

Endotoxins

Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test.

Tested

ELISA, MS, SDS-PAGE, WB

Applications
Target

NTNG1

Preservatives

Lyophilized from a 0.2 μM filtered solution of 20mM PB, 150mM NaCl, pH 7.4.

Form/Appearance

Lyophilized from a 0.2 μM filtered solution of 20mM PB, 150mM NaCl, pH 7.4.

Storage

Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 2-8°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.

Note

For research use only

Application notes

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 μg/ml. Dissolve the lyophilized protein in 1X PBS. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

Protein Sequence

HYDLCKTQIY TEEGKVWDYM ACQPESTDMT
KYLKVKLDPP DITCGDPPET FCAMGNPYMC
NNECDASTPE LAHPPELMFD FEGRHPSTFW
QSATWKEYPK PLQVNITLSW SKTIELTDNI
VITFESGRPD QMILEKSLDY GRTWQPYQYY
ATDCLDAFHM DPKSVKDLSQ HTVLEIICTE
EYSTGYTTNS KIIHFEIKDR FAFFAGPRLR
NMAISLYGQLD TTKKLRDFFT VTDLRIRLLR
PAVGEIFVDE LHLARYFYAI SDIKVRGRCK
CNLHATVCVY DNSKLTCECE HNTTGPDCGK
CKKNYQGRPW SPGSYLPIPK GTANTCIPSI
SSIGTNVCDN ELLHCQNGGT CHNNVRCLCP
AAYTGILCEK LRCEEAGSCG SVDHHHHHHH

Purity

Greater than 95% as determined by reducing SDS-PAGE.

Source

Human cells