

Product Datasheet

Human PACSIN2 protein (orb392230)

Description

Recombinant human PACSIN2 protein is

system and the target gene...

Reactivity

Human

Endotoxins

Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test.

Tested

ELISA, MS, SDS-PAGE, WB

Applications

Target

PACSIN2

Preservatives

Lyophilized from a 0.2 μM filtered solution of 20mM PB, 150mM NaCl, pH 7.4.

Form/Appearance

Lyophilized from a 0.2 μM filtered solution of 20mM PB, 150mM NaCl, pH 7.4.

Storage

Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 2-8°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.

Note

For research use only

Application notes

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 μg/ml. Dissolve the lyophilized protein in 1X PBS. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

Protein Sequence

MSVTYDDSVG VEVSSDSFWE VGNYKRTVKR
IDDGHRCLSD LMNCLHERAR IEKAYAQQLT
EWARRWRQLV EKG PQYGTVE KAWMAFMSEA
ERVSELHLEV KASLMNDDFE KIKNWQKEAF
HKQMMGGFKE TKEAEDGFRK AQPWAKKLK
EVEAAKKAHH AACKEEKLAI SREANSKADP
SLNPEQLKKL QDKIEKCKQD VLKTEKEYEK
SLKELDQGTP QYMENMEQVF EQCQQFEEKR
LRFFREVILLE VQKHLDSLNV AGYKAIYHDL
EQSIRAADAV EDLRWFRANH GPGMAMNWPQ
FEEWSADLNR TLRREKKKA TDGVTLTGIN
QTGDQSLPSK PSSTLNVP SN PAQSAQSQSS
YNPFEDEDDT GSTVSEKDDT KAKNVSSYEK
TQSYPTDWS D DESNNPFSST DANGDSNPFD
DDATSGTEVR VRALYDIEGQ EHDELSFKAG
DELTKMEDED EQGWCKGR LD NGQVGLYPAN
YVEAIQVDHH HHHH

Biorbyt Ltd.

7 Signet Court, Swann's Road, Cambridge, CB5 8LA, United Kingdom

Email: info@biorbyt.com | Phone: +44 (0) 1223 859-353 | Fax: +44 (0)1223 280 240

Biorbyt LLC.

68 TW Alexander Drive
Research Triangle Park
Durham, North Carolina
27709. United States

Email: info@biorbyt.com | Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558