



Product Datasheet

ACTR3B antibody (orb39227)

Description	Rabbit polyclonal antibody to ACTR3B
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	ACTR3B
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	47608
Uniprot ID	Q9P1U1
NCBI	NP_065178.1 , NP_001035225.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.



Western blot analysis of ACTR3B antibody...