

Product Datasheet

Human PIGR protein (orb392281)

Description

Recombinant human PIGR protein is produced

target gene en...

Reactivity

Human

Endotoxins

Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test.

Tested

ELISA, MS, SDS-PAGE, WB

Applications
Target

PIGR

Preservatives

Lyophilized from a 0.2 μM filtered solution of 20mM PB, 150mM NaCl, pH 7.2.

Form/Appearance

Lyophilized from a 0.2 μM filtered solution of 20mM PB, 150mM NaCl, pH 7.2.

Storage

Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 2-8°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.

Note

For research use only

Application notes

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 μg/ml. Dissolve the lyophilized protein in 1X PBS. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

Protein Sequence

KSPIFGPEEV NSVEGNSVSI TCYYPPTSVN
 RHTRKYWCRQ GARGGCITLI SSEGIVSSKY
 AGRANLTNFP ENGTFVFNIA QLSQDDSGRY
 KCGLGINSRG LSFDSLEVS QGPGLLNDTK
 VYTVDLGRTV TINCPFKTEN AQKRKSLYKQ
 IGLYPVLVID SSGYVNPNTY GRIRLDIQGT
 GQLLSVVIN QLRLSDAGQY LCQAGDDSNS
 NKKNADLQVL KPEPELVYED LRGSVTFHCA
 LGPEVANVAK FLCRQSSGEN CDVVVNTLGK
 RAPAFEGRIL LNPQDKDSF SVVITGLRKE
 DAGRYLCGAH SDGQLQEGSP IQAWQLFVNE
 ESTIPRSPTV VKGVAGGSVA VLCYPNRKES
 KSIKYWCLWE GAQNGRCPLL VDSEGWVKAQ
 YEGRLSLEE PGNGTFTVIL NQLTSRDAGF
 YWCLTNGDTL WRTTVEIKII EGEPNLKVPK
 NVTAVLGETL KVPCHFCKF SSYEKYWCKW
 NNTGCQALPS QDEGPSKAFV NCDENSRLVS
 LTI NI VTRAD FGWYWCYVKO GHEYGETAAV

Biorbyt Ltd.

7 Signet Court, Swann's Road, Cambridge, CB5 8LA, United Kingdom

Email: info@biorbyt.com | Phone: +44 (0) 1223 859-353 | Fax: +44 (0)1223 280 240

Biorbyt LLC.

68 TW Alexander Drive
Research Triangle Park
Durham, North Carolina
27709, United States

Email: info@biorbyt.com | Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558