

Product Datasheet

Human PITPNA protein (orb392287)

Description

Recombinant human PITPNA protein is

and the target gene en...

Reactivity

Human

Endotoxins

Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test.

Tested

ELISA, MS, SDS-PAGE, WB

Applications

Target

PITPNA

Preservatives

Liquid in a 0.2 μM filtered solution of 20mM Tris-HCl, 1mM EDTA, 1mM DTT, pH 8.0.

Form/Appearance

Liquid in a 0.2 μM filtered solution of 20mM Tris-HCl, 1mM EDTA, 1mM DTT, pH 8.0.

Storage

Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 μg/ml. Dissolve the lyophilized protein in 1X PBS. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

Protein Sequence

MGSSHHHHHH SSGLVPRGSH MVLLKEYRVI
LPVSVDEYQV GQLYSVAEAS KNETGGGEGV
EVLVNEPYEK DGEKGQYTHK IYHLQSKVPT
FVRMLAPEGA LNIHEKAWNA YPYCRTVITN
EYMKEDFLIK IETWHKPD LG TQENVHKLEP
EAWKHVEAVY IDIADRSQVL SKDYKAEEDP
AKFKSIKTGR GPLGPNWKQE LVNQKDCPYM
CAYKLTVKF KWWGLQNKVE NFIHKQERRL
FTNFHRQLFC WLDKWVDLTM DDIRRMEEET
KRQLDEMQRK DPVKGMTADD

Purity

Greater than 95% as determined by reducing SDS-PAGE.

Source

E. coli

Uniprot ID

[Q00169](#)

Entrez

[5306](#)

Expiration Date

6 months from date of receipt.