

Product Datasheet

Human PMM1 protein (orb392300)

Description

Recombinant human PMM1 protein is

and the target gene encod...

Reactivity

Human

Endotoxins

Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test.

Tested

Applications

ELISA, MS, SDS-PAGE, WB

Target

PMM1

Preservatives

Liquid in a 0.2 μM filtered solution of 20mM Tris-HCl, 150mM NaCl, 1mM DTT, pH 8.0.

Form/Appearance

Liquid in a 0.2 μM filtered solution of 20mM Tris-HCl, 150mM NaCl, 1mM DTT, pH 8.0.

Storage

Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 μg/ml. Dissolve the lyophilized protein in 1X PBS. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

Protein Sequence

MAVTAQAARR KERVLCFLDV DGTLPARQK
IDPEVAAFLQ KLRSRVQIGV VGGSDYCKIA
EQLGDGDEVI EKFDYVFAEN GTVQYKHGRL
LSKQTIQNHL GEELLQDLIN FCLSYMALLR
LPKKRGTFIE FRNGMLNISP IGRSCTLEER
IEFSELDKKE KIREKFVEAL KTEFAGKGLR
FSRGGMISFD VFPEGWDKRY CLDSLDQDSF
DTIHFFGNET SPGGNDFEIF ADPRTVGHSV
VSPQDTVQRC REIFFPETAH EAVEHHHHHH

Purity

Greater than 95% as determined by reducing SDS-PAGE.

Source

E. coli

Uniprot ID

[Q92871](#)

Entrez

[5372](#)

Expiration Date

6 months from date of receipt.