

Product Datasheet

Human PRCP protein (orb392328)

Description

Recombinant human PRCP protein is

system and the target gene en...

Reactivity

Human

Endotoxins

Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test.

Tested

ELISA, MS, SDS-PAGE, WB

Applications

Target

PRCP

Preservatives

Liquid in a 0.2 μM filtered solution of 20mM NaAc-HAc, 150mM NaCl, 10% Glycerol, pH 4.5.

Form/Appearance

Liquid in a 0.2 μM filtered solution of 20mM NaAc-HAc, 150mM NaCl, 10% Glycerol, pH 4.5.

Storage

Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 μg/ml. Dissolve the lyophilized protein in 1X PBS. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

Protein Sequence

LRPALRALGS LHLPTNPTSL PAVAKNYSVL
YFQQKVDHFG FNTVKTFNQR YLVADKYWKK
NGGSILFYTG NEGDIIWFCN NTGFMWDVAE
ELKAMLVFAE HRYYGESLPF GDNSFKDSRH
LNFLTSEQAL ADFAELIKHL KRTIPGAENQ
PVIAIGGSYG GMLAAWFRMK YPHMVVGALA
ASAPIWQFED LVPCGVFMKI VTDFRKSGP
HCSESIHRWS DAINRLSNTG SGLQWLTGAL
HLCSPLOTSQD IQHLKDWIS TWVNLAMVDY
PYASNFLQPL PAWPIKVVCQ YLKNPNVSDS
LLLQNIFQAL NVYYNYSQGV KCLNISETAT
SSLGTLGWSY QACTEVVMPF CTNGVDDMFE
PHSWNLKELS DDCFQQWGV PRPSWITMY
GGKNISHTN IVFSNGELDP WSGGGVTKDI
TDTLVAVTIS EGAHHLDLRT KNALDPMSVL
LARSLEVRHM KNWIRDFYDS AGKQHVDDHHH
HHH

Biorbyt Ltd.

7 Signet Court, Swann's Road, Cambridge, CB5 8LA, United Kingdom

Email: info@biorbyt.com | Phone: +44 (0) 1223 859-353 | Fax: +44 (0)1223 280 240

Biorbyt LLC.

68 TW Alexander Drive
Research Triangle Park
Durham, North Carolina
27709, United States

Email: info@biorbyt.com | Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558