

Product Datasheet

Human PSMD10 protein (orb392354)

Description

Recombinant human PSMD10 protein is

and the target gene en...

Reactivity

Human

Endotoxins

Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test.

Tested

ELISA, MS, SDS-PAGE, WB

Applications

Target

PSMD10

Preservatives

Liquid in a 0.2 μM filtered solution of PBS, 10% Glycerol, pH 7.4.

Form/Appearance

Liquid in a 0.2 μM filtered solution of PBS, 10% Glycerol, pH 7.4.

Storage

Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 μg/ml. Dissolve the lyophilized protein in 1X PBS. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

Protein Sequence

MEGCVSNLMV CNLAYSGKLE ELKESILADK
SLATRTDQDS RTALHWACSA GHTEIVEFLL
QLGVPVNDKD DAGWSPLHIA ASAGRDEIVK
ALLGKGAQVN AVNQNGCTPL HYAASKNRHE
IAVMLLEGGA NPDAKDHYEA TAMHRAAAKG
NLKMIHILLY YKASTNIQDT EGNTPLHLAC
DEERVEEAKL LVSQGASIYI ENKEEKTPLQ
VAKGGLGLIL KRMVEG

Purity

Greater than 95% as determined by SEC-HPLC and reducing SDS-PAGE.

Source

E. coli

Uniprot ID

[075832](#)

Entrez

[5716](#)

Expiration Date

6 months from date of receipt.