

Product Datasheet

Human PVRL1 protein (orb392367)

Description

Recombinant human PVRL1 protein is

system and the target gene e...

Reactivity

Human

Endotoxins

Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test.

Tested

Applications

ELISA, MS, SDS-PAGE, WB

Target

PVRL1

Preservatives

Lyophilized from a 0.2 μM filtered solution of 20mM PB, 150mM NaCl, pH 7.2.

Form/Appearance

Lyophilized from a 0.2 μM filtered solution of 20mM PB, 150mM NaCl, pH 7.2.

Storage

Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 2-8°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.

Note

For research use only

Application notes

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 μg/ml. Dissolve the lyophilized protein in 1X PBS. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

Protein Sequence

QVVQVNSMY GFIGTDVVLH CSFANPLPSV
KITQVTWQKS TNGSKQNVAI YNPSMGVSVL
APYRERVEFL RPSFTDGTIR LSRLELEDEG
VYICEFATFP TGNRESQLNL TVMAKPTNWI
EGTQAVLRAK KGQDDKVLVA TCTSANGKPP
SVVSWETRLK GEAQEYQIRN PNGTVTVISR
YRLVPSREAH QQSLACIVNY HMDRFKESLT
LNVQYEPEVT IEGFDGNWYL QRMDVKLTCK
ADANPPATEY HWTTLNGLSLP KGVEAQNRTL
FFKGPINYSL AGTYICEATN PIGTRSGQVE
VNITVDHHHH HH

Purity

Greater than 95% as determined by SEC-HPLC and reducing SDS-PAGE.

Source

Human cells

Uniprot ID

Q15223