

Product Datasheet

Human PVRL4 protein (orb392371)

Description

Recombinant human PVRL4 protein is

system and the target gene e...

Reactivity

Human

Endotoxins

Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test.

Tested

ELISA, MS, SDS-PAGE, WB

Applications
Target

PVRL4

Preservatives

Lyophilized from a 0.2 μM filtered solution of 20mM PB, 150mM NaCl, pH 7.4.

Form/Appearance

Lyophilized from a 0.2 μM filtered solution of 20mM PB, 150mM NaCl, pH 7.4.

Storage

Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 2-8°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.

Note

For research use only

Application notes

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 μg/ml. Dissolve the lyophilized protein in 1X PBS. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

Protein Sequence

GELETSDVVT VVLGQDAKLP CFYRGDSGEQ
VGQVAWARVD AGEQAQELAL LHSKYGLHVS
PAYEGRVEQP PPRNPLDGS VLLRNAVQAD
EGEYECRVST FPAGSFQARL RLRVLVPLP
SLNPGPALEE GQGLTLAASC TAEGSPAPSV
TWDTEVKGTT SSRSFKHSRS AAVTSEFHLV
PSRSMNGQPL TCVVSHPGLL QDQRITHILH
VSFLAEASVR GLEDQNLWHI GREGAMLKCL
SEGQPPPSYN WTRLDGPLPS GVRVDGDTLG
FPPLTTEHSG IYVCHVSNEF SSRDSQVTVD
VLDPQEDSGK QVDLVASAVV VDHHHHHH

Purity

Greater than 95% as determined by reducing SDS-PAGE.

Source

Human cells

Uniprot ID
Q96NY8