

Product Datasheet

Human SDCBP protein (orb392427)

Description

Recombinant human SDCBP protein is

and the target gene enc...

Reactivity

Human

Endotoxins

Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test.

Tested

ELISA, MS, SDS-PAGE, WB

Applications

Target

SDCBP

Preservatives

Liquid in a 0.2 μM filtered solution of 20mM Tris HCl, 100mM NaCl, 40% glycerol, pH 8.0.

Form/Appearance

Liquid in a 0.2 μM filtered solution of 20mM Tris HCl, 100mM NaCl, 40% glycerol, pH 8.0.

Storage

Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 μg/ml. Dissolve the lyophilized protein in 1X PBS. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

Protein Sequence

SLYPSLEDLK VDKVIQAQTA FSANPANPAI
LSEASAPIPH DGNLYPRLYP ELSQYMGLSL
NEEEIRANVA VVSGAPLQGG LVARPSSINY
MVAPVTGNDV GIRRAEIKQG IREVILCKDQ
DGKIGLRLKS IDNGIFVQLV QANSPASLVG
LRFGDQVLQI NGENCAGWSS DKAHKVLKQA
FGEKITMTIR DRPFERTITM HKDSTGHVGF
IFKNGKITSI VKDSSAARNG LLTEHNICEI
NGQNVIGLKD SQIADILSTS GTVVTITIMP
AFIFEHIKR MAPSIMKSLM DHTIPEVLEH
HHHHH

Purity

Greater than 95% as determined by reducing SDS-PAGE.

Source

E. coli

Uniprot ID

[O00560](#)

Entrez

[6386](#)