

Product Datasheet

Human SENP2 protein (orb392445)

Description

Recombinant human SENP2 protein is

and the target gene enc...

Reactivity

Human

Endotoxins

Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test.

Tested

Applications

ELISA, MS, SDS-PAGE, WB

Target

SENP2

Preservatives

Liquid in a 0.2 μM filtered solution of 50mM HEPES, 5% glycerol, pH 7.4.

Form/Appearance

Liquid in a 0.2 μM filtered solution of 50mM HEPES, 5% glycerol, pH 7.4.

Storage

Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 μg/ml. Dissolve the lyophilized protein in 1X PBS. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

Protein Sequence

MGHHHHHHGS LQDSEVNQEA KPEVKPEVKP
ETHINLKVSD GSSEIFFKIK KTTPLRRLME
AFKRQ GKEM DSLRFLYDGI RIQADQAPED
LDMEDNDIIE AHREQIGGDD LLELTEDMEK
EISNALGHGP QDEILSSAFK LRITRGDIQT
LKNYHWLNDE VINFYMNLLV ERNKKQGYP
LHVFSTFFYP KLKSGGYQAV KRWTGKVNLF
EQEILVPIH RKVHWSLVVI DLRKKCLKYL
DSMGQKGHRI CEILLQYLQD ESKTKRNSDL
NLEWTHHSM KPHEIPQQLN GSDCGMFTCK
YADYISRDKP ITFTQHQMPL FRKKMVWEIL
HQQLL

Purity

Greater than 95% as determined by reducing SDS-PAGE.

Source

E. coli

Uniprot ID

[Q9HC62](#)

Entrez

[59343](#)

Biorbyt Ltd.

7 Signet Court, Swann's Road, Cambridge, CB5 8LA, United Kingdom

Email: info@biorbyt.com | Phone: +44 (0) 1223 859-353 | Fax: +44 (0)1223 280 240

Biorbyt LLC.

68 TW Alexander Drive
Research Triangle Park
Durham, North Carolina
27709. United States

Email: info@biorbyt.com | Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558