

Product Datasheet

Human SENP8 protein (orb392447)

Description

Recombinant human SENP8 protein is

and the target gene enc...

Reactivity

Human

Endotoxins

Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test.

Tested

ELISA, MS, SDS-PAGE, WB

Applications

Target

SENP8

Preservatives

Liquid in a 0.2 μM filtered solution of 20mM PB, 150mM NaCl, pH 7.4.

Form/Appearance

Liquid in a 0.2 μM filtered solution of 20mM PB, 150mM NaCl, pH 7.4.

Storage

Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 μg/ml. Dissolve the lyophilized protein in 1X PBS. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

Protein Sequence

MGSSHHHHHH SGLVPRGSH MDPVLSYMD
SLLRQSDVSL LDPPSWLNDH IIGFAFEYFA
NSQFHDCSDH VSFISPEVTQ FIKCTSNPAE
IAMFLEPLDL PNKRVVFLAI NDNSNQAAGG
THWSLLVYLQ DKNSFFHYDS HRSNSVHAK
QVAEKLEAFL GRKGDKLAFLV EEKAPAQONS
YDCGMYVICN TEALCQNFFR QQTESLLQLL
TPAYITKKRG EWKDLIATLA KK

Purity

Greater than 95% as determined by SEC-HPLC and reducing SDS-PAGE.

Source

E. coli

Uniprot ID

[Q96LD8](#)

Entrez

[123228](#)

Expiration Date

6 months from date of receipt.