

Product Datasheet

Human SOAT2 protein (orb392504)

Description

Recombinant human SOAT2 protein is

and the target gene enc...

Reactivity

Human

Endotoxins

Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test.

Tested

ELISA, MS, SDS-PAGE, WB

Applications

Target

SOAT2

Preservatives

Liquid in a 0.2 μM filtered solution of 20mM PB, 150mM NaCl, pH 7.4.

Form/Appearance

Liquid in a 0.2 μM filtered solution of 20mM PB, 150mM NaCl, pH 7.4.

Storage

Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 μg/ml. Dissolve the lyophilized protein in 1X PBS. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

Protein Sequence

MSDKIIHLTD DSFDTDVLKA DGAILVDFWA
EWCGPCKMIA PILDEIADEY QGKLTVAKLN
IDQNPGTAPK YGIRGIPTLL LFKNGEVAAT
KVGALSKGQL KEFLDANLAG SGSGHMHMHHH
HHSSGLVPRG SGMKETAAAK FERQHMDSPD
LGTDDDDKAM ADIGSMNAGS DPVVIVSAAR
TIIGSFNGAL AAVPVQDLGS TVIKEVLKRA
TVAPEDVSEV IFGHVLAAGC GQNPVRQASV
GAGIPYSVPA WSCQMICGSG LKAVCLAVQS
IGIGDSSIVV AGGMENMSKA PHLAYLRTGV
KIGEMPLTDS ILCDGLTDAF HNCHMGITAE
NVAKKWQVSR EDQDKVAVLS QNRTENAQKA
GHFDKEIVPV LVSTRRLGIE VKTDEFPRHG
SNIEAMSKLK PYFLTDGTGT VTPANASGIN
DGAAAVLMK KSEADKRGLT PLARIVWSQ
VGVEPSIMGI GPIPAIKQAV TKAGWSLEDV
DIFEINEAFA AVSAAIVKEL GLNPEKVNIE
GGAIALGHPL GASGCRILVT LLHTLERMGR
SRGVAALCIG GGMGIAMCVQ RE

Biorbyt Ltd.

7 Signet Court, Swann's Road, Cambridge, CB5 8LA, United Kingdom

Email: info@biorbyt.com | Phone: +44 (0) 1223 859-353 | Fax: +44 (0)1223 280 240

Biorbyt LLC.

68 TW Alexander Drive
Research Triangle Park
Durham, North Carolina
27709. United States

Email: info@biorbyt.com | Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558