

Product Datasheet

Human SPINK1 protein (orb392511)

Description

Recombinant human SPINK1 protein is

system and the target gene ...

Reactivity

Human

Endotoxins

Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test.

Tested Applications

ELISA, MS, SDS-PAGE, WB

Target

SPINK1

Preservatives

Liquid in a 0.2 μM filtered solution of 20mM MES, 150mM NaCl, 2mM CaCl₂, 1mM DTT, 0.05% Brij35, 10% Glycerol, pH 6.0.

Form/Appearance

Liquid in a 0.2 μM filtered solution of 20mM MES, 150mM NaCl, 2mM CaCl₂, 1mM DTT, 0.05% Brij35, 10% Glycerol, pH 6.0.

Storage

Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 μg/ml. Dissolve the lyophilized protein in 1X PBS. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

Protein Sequence

DSLGREAKCY NELNGCTKIY DPVCGTDGNT YPNECVLCFE NRKRQTSILI QKSGPCVDHH HHHH

Purity

Greater than 95% as determined by SEC-HPLC and reducing SDS-PAGE.

Source

Human cells

Uniprot ID
P00995
Entrez
6690
Expiration Date

6 months from date of receipt.