

Product Datasheet

Mouse SPOCK3 protein (orb392515)

Description

Recombinant mouse SPOCK3 protein is

system and the target gene ...

Reactivity

Mouse

Endotoxins

Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test.

Tested Applications

ELISA, MS, SDS-PAGE, WB

Target

SPOCK3

Preservatives

Lyophilized from a 0.2 μM filtered solution of 20mM PB, 150mM NaCl, pH 7.4.

Form/Appearance

Lyophilized from a 0.2 μM filtered solution of 20mM PB, 150mM NaCl, pH 7.4.

Storage

Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 2-8°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.

Note

For research use only

Application notes

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 μg/ml. Dissolve the lyophilized protein in 1X PBS. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

Protein Sequence

AAAVAVAGGR SDGGNFLDEK QWLTTSQYD
KEVGQWNKFR DEVEDDYFRT WNP GKPF DQA
LDP AKDPCLK TKCSRHK VCI TQDAQTALCI
SHRRLTHSMK EVGGSHKQWR GLPSSTCKPC
PIAYASPVC G SDGHSYSSQC KLEYQACVLG
KQISIKCEGR CPCPSDKSMN IGRNVKRACS
DLEFREVANR LRDWFKALHE SGSQNKKTKA
LLRPERSRFD TSILPICKDS LGWMFNRLDT
NYDLLLDQSE LGSYLDKNE QCTKAFFNSC
DTYKDSLISN NEWCYCFQRQ QDPPCHTELS
NIQKRQGIKK LLGQYIPLCD EDGYKPTQC
HGSVGQCWCV DRYGNEVVG S RINGVADCAI
DFEISGDFAS GDFREWTDDE GEEDDIMNDK
DDIEDDDEDE GDDDDGDVH DGYIVDHHHH
HH

Purity

Greater than 95% as determined by reducing SDS-PAGE.

Biorbyt Ltd.

7 Signet Court, Swann's Road, Cambridge, CB5 8LA, United Kingdom

Email: info@biorbyt.com | Phone: +44 (0) 1223 859-353 | Fax: +44 (0)1223 280 240

Biorbyt LLC.

68 TW Alexander Drive
Research Triangle Park
Durham, North Carolina
27709. United States

Email: info@biorbyt.com | Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558

Expiration Date

6 months from date of receipt.