

Product Datasheet

Human STX8 protein (orb392528)

Description

Recombinant human STX8 protein is

and the target gene enco...

Reactivity

Human

Endotoxins

Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test.

Tested

Applications

ELISA, MS, SDS-PAGE, WB

Target

STX8

Preservatives

Lyophilized from a 0.2 μM filtered solution of 20mM PB, 150mM NaCl, pH 7.4.

Form/Appearance

Lyophilized from a 0.2 μM filtered solution of 20mM PB, 150mM NaCl, pH 7.4.

Storage

Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 2-8°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.

Note

For research use only

Application notes

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 μg/ml. Dissolve the lyophilized protein in 1X PBS. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

Protein Sequence

MAPDPWFSTY DSTCQIAQEI AEKIQQRNQY
ERKGEKAPKL TVTIRALLQN LKEKIALLLKD
LLLRAVSTHQ ITQLEGDRRQ NLLDDLVTRE
RLLASFKNE GAEPDLIRSS LMSEEAKRGA
PNPWLFEPE ETRGLGFDEI RQQQKIIQE
QDAGLDALSS IISRQKQMGQ EIGNELDEQN
EIIDDLANLV ENTDEKLRNE TRRVNMVDRK
SASCG

Purity

Greater than 95% as determined by reducing SDS-PAGE.

Source

E. coli

Uniprot ID

[Q9UNK0](#)

Entrez

[9482](#)

Expiration Date

6 months from date of receipt.