

Product Datasheet

Human THOP1 protein (orb392561)

Description

Recombinant human THOP1 protein is

and the target gene enc...

Reactivity

Human

Endotoxins

Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test.

Tested

ELISA, MS, SDS-PAGE, WB

Applications
Target

THOP1

Preservatives

Liquid in a 0.2 μM filtered solution of 20mM PB, 100mM NaCl, 10% Glycerol, pH 7.0.

Form/Appearance

Liquid in a 0.2 μM filtered solution of 20mM PB, 100mM NaCl, 10% Glycerol, pH 7.0.

Storage

Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 μg/ml. Dissolve the lyophilized protein in 1X PBS. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

Protein Sequence

MKPPAACAGD MADAASPCSV VNDLRWDLSA
 QQIEERTREL IEQTKRVDYD VGTQEFEDVS
 YESTLKALAD VEVYTVQRN ILDFPQHVS
 SKDIRASTE ADKKLSEFDV EMSMREDVYQ
 RIVWLQEKVQ KDSLRPEAAR YLERLIKLR
 RNGLHLPRET QENIKRIKKK LSLCIDFNK
 NLNEDTTFLP FTLQELGGLP EDFLNSLEKM
 EDGKLKVTLK YPHYFLLKK CHVPETRRKV
 EEFNCRCKE ENCAILKELV TLRAQKSRL
 GFHTHADYVL EMNMAKTSQT VATFLDELAQ
 KKLPLGEQER AVILELKRAE CERRGLPFDG
 RIRAWDMRYM MNQVEETRYC VDQNLKEYF
 PVQVVTHGLL GIYQELLGLA FHHEEGASAW
 HEDVRLYTAR DAASGEVVGK FYLDLYPREG
 KYGHAACFGL QPGCLRQDGS RQIAIAAMVA
 NFKPTADAP SLLQHDEVET YFHEFGHVMH
 QLCSQAEFAM FSGTHVERDF VEAPSQMLEN
 WVWEQEPLLR MSRHYRTGSA VPRELLEKLI
 ESRQANTGLF NLRQIVLAKV DQALHTQTD
 DPAFFYARLC OFILGVPATP GTNMPATEGH

Biorbyt Ltd.

7 Signet Court, Swann's Road, Cambridge, CB5 8LA, United Kingdom

Email: info@biorbyt.com | Phone: +44 (0) 1223 859-353 | Fax: +44 (0)1223 280 240

Biorbyt LLC.

68 TW Alexander Drive
Research Triangle Park
Durham, North Carolina
27709. United States

Email: info@biorbyt.com | Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558