

Product Datasheet

Mouse TNFRSF10B protein (orb392585)

Description

Recombinant mouse TNFRSF10B protein is

system and the target ge...

Reactivity

Mouse

Endotoxins

Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test.

Tested Applications

ELISA, MS, SDS-PAGE, WB

Target

TNFRSF10B

Preservatives

Lyophilized from a 0.2 μM filtered solution of 20mM PB, 150mM NaCl, pH 7.4.

Form/Appearance

Lyophilized from a 0.2 μM filtered solution of 20mM PB, 150mM NaCl, pH 7.4.

Storage

Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 2-8°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.

Note

For research use only

Application notes

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 μg/ml. Dissolve the lyophilized protein in 1X PBS. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

Protein Sequence

NPAHNRPAGL QRPEESPSRG PCLAGQYLSE
GNCKPCREGI DYTSHSNHSL DSCILCTVCK
EDKVVETRCN ITTNTVCRCK PGT FEDKDSP
EICQSCSNCT DGEEELTSTCT PRENRKCVSK
TAWASVDDIE GRMDEPKSCD KTHTCPPCPA
PELLGGPSVF LFPPKPKDTL MISRTPEVTC
VVVDVSHEDP EVKFNWYVDG VEVHNAKTKP
REEQYNSTYR VVSVLTVLHQ DWLNGKEYKC
KVSNAKALPAP IEKTISKAKG QPREPQVYTL
PPSREEMTKN QVSLTCLVKG FYPSDIAVEW
ESNGQPENNY KTTTPVLDSD GSFFLYSKLT
VDKSRWQQGN VFSCSVMHEA LHNHYTQKSL
SLSPGK

Purity

Greater than 95% as determined by reducing SDS-PAGE.

Source

Human cells