

Product Datasheet

Mouse TNFRSF12A protein (orb392591)

Description

Recombinant mouse TNFRSF12A protein is

system and the target ge...

Reactivity

Mouse

Endotoxins

Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test.

Tested Applications

ELISA, MS, SDS-PAGE, WB

Target

TNFRSF12A

Preservatives

Lyophilized from a 0.2 μM filtered solution of PBS, pH 7.4.

Form/Appearance

Lyophilized from a 0.2 μM filtered solution of PBS, pH 7.4.

Storage

Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 2-8°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.

Note

For research use only

Application notes

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 μg/ml. Dissolve the lyophilized protein in 1X PBS. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

Protein Sequence

EQAPGTSPCS SGSSWSADLD KCMDCASCPA
RPHSDFCLGC AAAPPAHFRL LWVDDIEGRM
DEPKSCDKTH TCPPCPAPEL LGGPSVFLFP
PKPKDTLMIS RTPEVTCVVV DVSHEDPEVK
FNWYVDGVEV HNAKTKPREE QYNSTYRVVS
VLTVLHQDWL NGKEYKCKVS NKALPAPIEK
TISKAKGQPR EPQVYTLPPS REEMTKNQVS
LTCLVKGFYP SDIAVEWESN GQPENNYKTT
PPVLDSGGSF FLYSKLTVDK SRWQQGNVFS
CSVMHEALHN HYTKSLSLS PGK

Purity

Greater than 95% as determined by reducing SDS-PAGE.

Source

Human cells

Uniprot ID

[Q9CR75](#)

Entrez

[27279](#)