

Product Datasheet

Mouse TNFRSF13B protein (orb392592)

Description

Recombinant mouse TNFRSF13B protein is

system and the target ge...

Reactivity

Mouse

Endotoxins

Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test.

Tested

Applications

ELISA, MS, SDS-PAGE, WB

Target

TNFRSF13B

Preservatives

Lyophilized from a 0.2 μM filtered solution of PBS, pH 7.4.

Form/Appearance

Lyophilized from a 0.2 μM filtered solution of PBS, pH 7.4.

Storage

Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 2-8°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.

Note

For research use only

Application notes

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 μg/ml. Dissolve the lyophilized protein in 1X PBS. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

Protein Sequence

MFCPKDQYWD SSRKSCVSCA LTCSQRSQRT
CTDFCKFINC RKEQGRIYDH LLGACVSCDS
TCTQHPQCA HFCEKRPRSQ ANLQPELGRP
QAGEVEVRSD NSGRHQGSEH GPGLRLSSDQ
LTLYCTVDDI EGRMDEPKSC DKTHTCPPCP
APELLGGPSV FLFPPKPKDT LMISRTPEVT
CVVVDVSHED PEVKFNWYVD GVEVHNAKTK
PREEQYNSTY RVVSVLTVLH QDWLNGKEYK
CKVSNKALPA PIEKTISKAK GQPREPQVYT
LPPSREEMTK NQVSLTCLVK GFYPSDIAVE
WESNGQPENN YKTTTPVLDS DGSFFLYSKL
TVDKSRWQQG NVFSCSVMHE ALHNHYTQKS
LSLSPGK

Purity

Greater than 95% as determined by reducing SDS-PAGE.

Source

Human cells

Biorbyt Ltd.

7 Signet Court, Swann's Road, Cambridge, CB5 8LA, United Kingdom

Email: info@biorbyt.com | Phone: +44 (0) 1223 859-353 | Fax: +44 (0)1223 280 240

Biorbyt LLC.

68 TW Alexander Drive
Research Triangle Park
Durham, North Carolina
27709, United States

Email: info@biorbyt.com | Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558