

Product Datasheet

Mouse TNFRSF1B protein (orb392600)

Description

Recombinant mouse TNFRSF1B protein is

system and the target gen...

Reactivity

Mouse

Endotoxins

Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test.

Tested

Applications

ELISA, MS, SDS-PAGE, WB

Target

TNFRSF1B

Preservatives

Lyophilized from a 0.2 μM filtered solution of PBS, pH 7.4.

Form/Appearance

Lyophilized from a 0.2 μM filtered solution of PBS, pH 7.4.

Storage

Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 2-8°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.

Note

For research use only

Application notes

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 μg/ml. Dissolve the lyophilized protein in 1X PBS. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

Protein Sequence

VPAQVVLTPY KPEPGYECQI SQEYYDRKAQ
MCCAKCPPGQ YVKHFCNKTS DTVCADCEAS
MYTQVWNQFR TCLSCSSSCT TDQVEIRACT
KQQNRVCACE AGRYCALKTH SGSCRQCMRL
SKCGPGFGVA SSRAPNGNVL CKACAPGTFS
DTTSSTDVCR PHRICSILAI PGNASTDAVC
APESPTLSAI PRTLVSQPE PTRSQPLDQE
PGPSQTPSIL TSLGSTPIIE QSTKGGVDHH
HHHH

Purity

Greater than 95% as determined by reducing SDS-PAGE.

Source

Human cells

Uniprot ID

P25119

Entrez

21938