

Product Datasheet

Human TNFRSF25 protein (orb392604)

Description

Recombinant human TNFRSF25 protein is

system and the target gen...

Reactivity

Human

Endotoxins

Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test.

Tested

Applications

ELISA, MS, SDS-PAGE, WB

Target

TNFRSF25

Preservatives

Lyophilized from a 0.2 μM filtered solution of PBS, pH 7.4.

Form/Appearance

Lyophilized from a 0.2 μM filtered solution of PBS, pH 7.4.

Storage

Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 2-8°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.

Note

For research use only

Application notes

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 μg/ml. Dissolve the lyophilized protein in 1X PBS. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

Protein Sequence

QGGTRSPRCD CAGDFHKKIG LFCCRGCPAG
HYLKAPCTEP CGNSTCLVCP QDTFLAWENH
HNSECARCQA CDEQASQVAL ENCSAVADTR
CGCKPGWFVE CQVSQCVSSS PFYCQPCLDC
GALHRHTRLL CSRRDTCGT CLPGFYEHD
GCVSCPTSTL GSCPERCAAV CGWRQMFVDD
IEGRMDEPKS CDKTHTCPPC PAPELLGGPS
VFLFPPKPKD TLMISRTPEV TCVVVDVSHE
DPEVKFNWYV DGVEVHNAKT KPREEQYNST
YRVVSVLTVL HQDWLNGKEY KCKVSNKALP
APIEKTISKA KGQPREPQVY TLPPSREEMT
KNQVSLTCLV KGFYPSDIAV EWESNGQPEN
NYKTTTPVLD SDGSFFLYSK LTVDKSRWQQ
GNVFSCSVMH EALHNHYTQK SLSLSPGK

Purity

Greater than 95% as determined by reducing SDS-PAGE.