

Product Datasheet

Human TNFRSF4 protein (orb392606)

Description

Recombinant human TNFRSF4 protein is

system and the target gene...

Reactivity

Human

Endotoxins

Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test.

Tested

Applications

ELISA, MS, SDS-PAGE, WB

Target

TNFRSF4

Preservatives

Lyophilized from a 0.2 μM filtered solution of PBS, pH 7.4.

Form/Appearance

Lyophilized from a 0.2 μM filtered solution of PBS, pH 7.4.

Storage

Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 2-8°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.

Note

For research use only

Application notes

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 μg/ml. Dissolve the lyophilized protein in 1X PBS. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

Protein Sequence

DKTHTCPPCP APELLGGPSV FLFPPKPKDT
LMISRTPEVT CVVVDVSHED PEVKFNWYVD
GVEVHNAKTK PREEQYNSTY RVVSVLTVLH
QDWLNGKEYK CKVSNKALPA PIEKTISKAK
GQPREPQVYT LPPSRDELTK NQVSLTCLVK
GFYPSDIAVE WESNGQPENN YKTTTPVLDS
DGSFFLYSKL TVDKSRWQQG NVFSCSV MHE
ALHNHYTQKS LSLSPGKDDD DKLHC VGDY
PSNDRCCHEC RPGNGMVSRC SRSQNTVCRP
CGPGFYNDVV SSKPCKPCTW CNLRSGSERK
QLCTATQDTV CRCRAGTQPL DSYKPGVDCA
PCPPGHFSPG DNQACKPWTN CTLAGKHTLQ
PASNSSDAIC EDRDPPATQP QETQGPPARP
ITVQPTEAWP RTSQGPSTRP VEVPGGRAVA

Purity

Greater than 95% as determined by reducing SDS-PAGE.