

Product Datasheet

Human TNFRSF9 protein (orb392611)

Description

Recombinant human TNFRSF9 protein is

system and the target gene...

Reactivity

Human

Endotoxins

Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test.

Tested Applications

ELISA, MS, SDS-PAGE, WB

Target

TNFRSF9

Preservatives

Lyophilized from a 0.2 μM filtered solution of PBS pH 7.4.

Form/Appearance

Lyophilized from a 0.2 μM filtered solution of PBS pH 7.4.

Storage

Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 2-8°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.

Note

For research use only

Application notes

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 μg/ml. Dissolve the lyophilized protein in 1X PBS. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

Protein Sequence

LQDPCSNCPA GTFCDNNRNQ ICSPCPPNSF
SSAGGQRTCD ICRQCKGVFR TRKECSSTSN
AECDCPTGFH CLGAGCSMCE QDCKQGQELT
KKGCKDCCFG TFNDQKRGIC RPWTNCSLDG
KSVLVNGTKE RDVVCGPSA DLSPGASSVT
PPAPAREPGH SPQGGGSDK THTCPPCPAP
ELLGGPSVFL FPPKPKDTLM ISRTPEVTCV
VVDVSHEDPE VKFNWYVDGV EVHNAKTKPR
EEQYNSTYRV VSVLTVLHQD WLNGKEYKCK
VSNKALPAPI EKTISKAKGQ PREPQVYTL
PSREEMTKNQ VSLTCLVKGF YPSDIAVEWE
SNGQPENNYK TTPPVLDSDG SFFLYSKLTV
DKSRWQQGNV FSCVMHEAL HNHYTQKSLS
LSPGKHHHHH H

Purity

Greater than 95% as determined by reducing SDS-PAGE.