

Product Datasheet

Mouse TNFRSF9 protein (orb392612)

Description

Recombinant mouse TNFRSF9 protein is

system and the target gene...

Reactivity

Mouse

Endotoxins

Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test.

Tested

Applications

ELISA, MS, SDS-PAGE, WB

Target

TNFRSF9

Preservatives

Liquid in a 0.2 μM filtered solution of 20mM PB, 150mM NaCl, pH 7.4.

Form/Appearance

Liquid in a 0.2 μM filtered solution of 20mM PB, 150mM NaCl, pH 7.4.

Storage

Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 μg/ml. Dissolve the lyophilized protein in 1X PBS. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

Protein Sequence

VQNSCDNCQP GTFCRKYNPV CKSCPPSTFS
SIGGQPNCNI CRVCAGYFRF KKFCSSSTHNA
ECECIEGFHC LGPQCTRCEK DCRPGQELTK
QGCKTCSLGT FNDQNGTGVC RPWTNCSLDG
RSVLKTGTTE KDVVCGPPVV SFSPSTTISV
TPEGPGGHS LQVLVDDIEG RMDEPKSCDK
THTCPPCPAP ELLGGPSVFL FPPKPKDTLM
ISRTPEVTCV VVDVSHEDPE VKFNWYVDGV
EVHNAKTKPR EEQYNSTYRV VSVLTVLHQD
WLNQKEYKCK VSNKALPAPI EKTISKAKGQ
PREPQVYTLPSREEMTKNQ VSLTCLVKGF
YPSDIAVEWE SNGQPENNYK TTPPVLDSDG
SFFLYSKLTV DKSRWQQGNV FSCSVMEAL
HNHYTQKSLS LSPGK

Purity

Greater than 95% as determined by reducing SDS-PAGE.

Source

Human cells

Biorbyt Ltd.

7 Signet Court, Swann's Road, Cambridge, CB5 8LA, United Kingdom

Email: info@biorbyt.com | Phone: +44 (0) 1223 859-353 | Fax: +44 (0)1223 280 240

Biorbyt LLC.

68 TW Alexander Drive
Research Triangle Park
Durham, North Carolina
27709, United States

Email: info@biorbyt.com | Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558