

Product Datasheet

Human TXNDC12 protein (orb392645)

Description

Recombinant human TXNDC12 protein is

system and the target gene...

Reactivity

Human

Endotoxins

Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test.

Tested

ELISA, MS, SDS-PAGE, WB

Applications

Target

TXNDC12

Preservatives

Liquid in a 0.2 μM filtered solution of 20mM PB, 150mM NaCl, 10% Glycerol, pH 7.4.

Form/Appearance

Liquid in a 0.2 μM filtered solution of 20mM PB, 150mM NaCl, 10% Glycerol, pH 7.4.

Storage

Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 μg/ml. Dissolve the lyophilized protein in 1X PBS. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

Protein Sequence

HNGLGKGFGD HIHWRTLEDG KKEAAASGLP
LMVIIHKS WC GACKALKPKF AESTEISELS
HNFVMVNLED EEEKDEDFD PDGGYIPRIL
FLDPGKGVHP EIINENG NPS YKYFYVSAEQ
VVQGMKEAQE RLTGDAFRKK HLV DHHHHHHH

Purity

Greater than 95% as determined by reducing SDS-PAGE.

Source

Human cells

Uniprot ID

[O95881](#)

Entrez

[51060](#)

Expiration Date

6 months from date of receipt.