

Product Datasheet

Human TXNDC15 protein (orb392646)

Description

Recombinant human TXNDC15 protein is

system and the target gene...

Reactivity

Human

Endotoxins

Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test.

Tested

Applications

ELISA, MS, SDS-PAGE, WB

Target

TXNDC15

Preservatives

Lyophilized from a 0.2 μM filtered solution of 20mM PB, 150mM NaCl, pH 7.4.

Form/Appearance

Lyophilized from a 0.2 μM filtered solution of 20mM PB, 150mM NaCl, pH 7.4.

Storage

Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 2-8°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.

Note

For research use only

Application notes

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 μg/ml. Dissolve the lyophilized protein in 1X PBS. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

Protein Sequence

VEVAEESGRL WSEEQPAHPL QVGAVYLGEE
ELLHDPMGQD RAAEEANAVL GLDTQGDHNV
MLSVIPGEAE DKVSSEPSGV TCGAGGAEDS
RCNVRESLFS LDGAGAHFPD REEEYYTEPE
VAESDAAPTE DSNNTESLKS PKVNCEERNI
TGLENFTLKI LNMSQDLMDL LNPNGSDCTL
VLFYTPWCRF SASLAPHFNS LPRAFPAHF
LALDASQHSS LSTRFGTVAV PNILLFQGAK
PMARFNHTDR TLETLKIFIF NQTGIEAKKN
VVVTQADQIG PLPSTLIKSV DHHHHHH

Purity

Greater than 95% as determined by reducing SDS-PAGE.

Source

Human cells

Uniprot ID

[Q96J42](#)

Entrez

[79770](#)