

Product Datasheet

Mouse TYRO3 protein (orb392647)

Description

Recombinant mouse TYRO3 protein is

system and the target gene e...

Reactivity

Mouse

Endotoxins

Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test.

Tested

Applications

ELISA, MS, SDS-PAGE, WB

Target

TYRO3

Preservatives

Lyophilized from a 0.2 μM filtered solution of PBS, pH 7.4.

Form/Appearance

Lyophilized from a 0.2 μM filtered solution of PBS, pH 7.4.

Storage

Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 2-8°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.

Note

For research use only

Application notes

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 μg/ml. Dissolve the lyophilized protein in 1X PBS. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

Protein Sequence

AGLKLMGAPV KMTVSQGPV KLNCSVEGME
DPDIHWMKDG TVVQNASQVS ISISEHSWIG
LLSLKSVERS DAGLYWCQVK DGEETKISQS
VWLTVEGVPP FTVEPKDLAV PPNAPFQLSC
EAVGPPEPVT IYWWRGTLTKV GGPAPSPSVL
NVTGVTQRTE FSCEARNIKG LATSRAIVR
LQAPPAAPFN TTVTTISSYN ASVAWVPGAD
GLALLHSCTV QVAHAPGEWE ALAVVVPVPP
FTCLLRNLAP ATNYSRLVRC ANALGPSPYG
DWVPFQTKGL APARAPQNFH AIRTDSGLIL
EWEEVIPEDP GEGPLGPYKL SWVQENGTDQ
ELMVEGTRAN LTDWDPQKDL ILRVCASNAI
GDGPWSQPLV VSSHDHAGRQ GPPHSRTSGG
GGSGGGGSGG GGSPRDCGCK PCICTVPEVS
SVFIFPPKPK DVLITLTPK VTCVVVDISK
DDPEVQFSWF VDDVEVHTAQ TQPREEQFNS
TFRSVSELPI MHQDWLNGKE FKCRVNSAAF
PAPIEKTISK TKGRPKAPOV YTIPPPKFOM

Biorbyt Ltd.

7 Signet Court, Swann's Road, Cambridge, CB5 8LA, United Kingdom

Email: info@biorbyt.com | Phone: +44 (0) 1223 859-353 | Fax: +44 (0)1223 280 240

Biorbyt LLC.

68 TW Alexander Drive
Research Triangle Park
Durham, North Carolina
27709, United States

Email: info@biorbyt.com | Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558