

Product Datasheet

Human UBAP1 protein (orb392649)

Description

Recombinant human UBAP1 protein is

and the target gene enc...

Reactivity

Human

Endotoxins

Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test.

Tested

ELISA, MS, SDS-PAGE, WB

Applications

Target

UBAP1

Preservatives

Lyophilized from a 0.2 μM filtered solution of 20mM PB, 150mM NaCl, 1mM DTT, pH 7.4.

Form/Appearance

Lyophilized from a 0.2 μM filtered solution of 20mM PB, 150mM NaCl, 1mM DTT, pH 7.4.

Storage

Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 2-8°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.

Note

For research use only

Application notes

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 μg/ml. Dissolve the lyophilized protein in 1X PBS. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

Protein Sequence

MASKKLGADF HGTFSYLDDV PFKTGDKFKT
PAKVGLPIGF SLPDCLQVVR EVQYDFSLEK
KTIEWAEEIK KIEEAEREAE CKIAEAEAKV
NSKSGPEGDS KMSFSKTHST ATMPPPINPI
LASLQHNSIL TPTRVSSSAT KQKVLSPPHI
KADFNLADFE CEEDPFDNLE LKTIDEKEEL
RNILVGTTGP IMAQLLDNNL PRGGSGSVLQ
DEEVLASLER ATLDFKPLHK PNGFITLPQL
GNCEKMSLSS KVSLLPIPAV SNIKSLSPFK
LSDSDSNQKT AKLASTFHST SCLRNGTFQN
SLKPSTQSSA SELNGHHTLG LSALNLDST
EMPALTSSQM PSLSVLSVCT EESSPPNTGP
TVTPPNFSVS QVPNMPSCPQ AYSELQMLSP
SERQCVETVV NMGYSYECVL RAMKKKGNI
EQILDYLF AH GQLCEKGFDP LLVEEALEMH
QCSEKMMEF LQLMSKFEM GFELKDIKEV
LLHNNDQDN ALEDLMARAG ASLEHHHHHH

Biorbyt Ltd.

7 Signet Court, Swann's Road, Cambridge, CB5 8LA, United Kingdom

Email: info@biorbyt.com | Phone: +44 (0) 1223 859-353 | Fax: +44 (0)1223 280 240

Biorbyt LLC.

68 TW Alexander Drive
Research Triangle Park
Durham, North Carolina
27709, United States

Email: info@biorbyt.com | Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558