




Product Datasheet

alpha 2 Macroglobulin antibody (orb39265)

Description	Rabbit polyclonal antibody to alpha 2 Macroglobulin
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IHC-P, WB
Immunogen	Synthetic Peptide
Target	C5
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	188305
Uniprot ID	P01031
NCBI	NP_001726.2
Dilution Range	WB: 1:1000, IHC-P: 1:50-100, FACS: 1:10-50
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of formali...



Flow cytometric analysis of C5 antibody ...